



NANORIP

User Manual

For PRINTRONIX LP654C or LP844C printers

IMPORTANT:

This manual describes all the functionalities available in the full-option version of NANORIP Ver 1.0.0.x. Available functionalities of your product may be different from those described in this manual.

License, manual and trademark information :

RB3L and RB3L Logo are registered trademark of RB3L Ltd. RB3L Ltd owns all copyrights pertaining to this manual. Unauthorized reproduction, transcription, and translation of this document are prohibited. You must obtain written permission from RB3L Ltd to reproduce, transcribe, or translate any part of this manual. Copyright 2022 RB3L Ltd, All rights reserved.

NANORIP, Printing Program License Agreement.

The details of the software license agreement described below must be agreed by you before using the software.

Software License Agreement

IMPORTANT

Please read this Software License Agreement (Agreement) carefully before using RB3L NANORIP (Licensed Software) with a printer (Product).

This Licensed Software is intended to work only in conjunction with the Product

By using the Licensed Software with the Product, you are indicating your acceptance of these terms and conditions and this Agreement is deemed effective between you (Licensee) and RB3L Ltd (Licensor).

If you do not agree to the terms of this Agreement, you may not use the Licensed Software.

1. Scope of the license

Licensor grants, and Licensee accepts, a non-exclusive license to use the Licensed Software solely in conjunction with the Product.

2. Ownership and Restriction

(1) Licensor or its suppliers own all the copyright and proprietary rights in and to the Licensed Software. The structure, organization and code contained in Licensed Software are the valuable trade secrets of Licensor and its suppliers.

Licensed Software is also protected by United States Copyright Law and International Treaty provisions. Licensee must treat the Licensed Software just as Licensee would treat any other copyrighted material, such as a book.

(2) Licensee may not make copies rent, lease, distribute, transfer or reprint the Licensed Software, in whole or in part.

(3) Licensee agrees not to modify, alter, translate, reverse engineer, decompile, disassemble, extract in part or separate in part the Licensed Software.

(4) Licensee agrees not to change the file names for the Licensed Software.

(5) Except as stated in this Agreement, Licensor does not grant Licensee any intellectual property rights in or to the Licensed Software.

3. Term and Termination

(1) This Agreement is effective until the Product is destroyed.

(2) Licensor may terminate this Agreement if Licensee fails to comply with any of the terms and conditions of this Agreement. Upon termination, Licensee shall destroy Licensed Software and its copies in Licensee's possession and control.

4. Warranty

THE LICENSED SOFTWARE IS PROVIDED "AS IS". NEITHER LICENSOR NOR ITS SUPPLIERS WARRANT THAT THE OPERATION OF THE LICENSED SOFTWARE WILL BE UNINTERRUPTED, ERROR FREE, OR WILL MEET LICENSEES NEEDS. LICENSOR AND ITS SUPPLIERS MAKE NO WARRANTY, EXPRESS OR IMPLIED, AS TO NONINFRINGEMENT OF ANY THIRD PARTY'S RIGHTS, MERCHANTABILITY, OR FITNESS FOR ANY PARTICULAR PURPOSE.

5. Disclaimer of Liability

IN NO EVENT WILL LICENSOR OR ITS SUPPLIERS BE LIABLE TO LICENSEE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, CONSEQUENTIAL OR PUNITIVE DAMAGES, INCLUDING, BUT NOT LIMITED TO, LOSS OF PROFITS, LOSS OF SAVINGS, OR LOSS OF DATA, EVEN IF LICENSOR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGE NOR FOR ANY CLAIM BY ANY THIRD PARTY (IES) IN ANY WAY ARISING OUT OF OR RELATING TO THE LICENSED SOFTWARE, WHETHER THE CLAIM ALLEGES TORTIOUS CONDUCT (INCLUDING, BUT NOT LIMITED TO, NEGLIGENCE) OR ANY OTHER LEGAL THEORY, OR FOR ANY CLAIM BY ANY THIRD PARTY.

6. Governing Law

The license with respect to the Licensed Software will be governed by the laws in force in United Kingdom.

7. Severability

If any part of this Agreement is found void and unenforceable, it will not affect the validity of the balance of the Agreement, which shall remain valid and enforceable according to its terms.

8. Export Restriction

Licensee agrees that the Licensed Software will not be shipped, transferred, exported or re-exported into any country or used in any manner prohibited by U.S., Europe or any other applicable export laws or regulations. Licensee agrees that it will not export or re-export the Licensed Software or products produced therefrom in any form without appropriate U.S., Europe or any other governmental licenses.

This Agreement shall automatically terminate upon failure by Licensee to comply with this Section 8.

9. Entire Agreement

Licensee represents that Licensee reads and understands this Agreement and that this Agreement constitutes the entire understanding and agreement between Licensor and Licensee as to the license relating to the Licensed Software, and supersedes and replaces any prior agreement, written or oral.

The Licensees obligations in this Agreement constitute the obligations to Licensor and all owners of the right licensed to Licensee under this Agreement.

10. Notice to U.S. Government End Users

All Software provided to the U.S. Government pursuant to solicitations issued on or after December 1, 1995 is provided with the commercial license rights and restrictions described elsewhere herein. All Software provided to the U.S. Government pursuant to solicitations issued prior to December 1, 1995 is provided with "Restricted Rights" as provided for in FAR, 48 CFR 52.227-14 (JUNE 1987) or DFAR, 48 CFR 252.227-7013 (OCT 1988), as applicable.

"Software" in this section is deemed to be the Licensed Software defined in this Agreement.

Some of the Licensed Software may include a separate software license agreement and if you agree to the separate software license agreement, the terms of such agreement shall prevail for the use of the software.

11. Licensee acknowledges and agrees that Licensor's suppliers are third party beneficiaries of this Agreement, with the right to enforce the obligations set forth herein with respect to the respective technology of such suppliers and/or Licensor.

1. **Presentation :**

1.1 **Introduction**

NANORIP transforms PRINTRONIX LP654C or LP844C printers into DICOM imagers in order to produce high quality low-cost prints, on paper with different sizes (A6, A5, A4, A3, etc. ..) for all medical imaging modalities. This imager can generate both GSDF compatible grayscale or color pages, and match them to the screen rendering.

Fonctions

- * High print quality with ProQ2400© multi-level technology
- * DICOM 3.0 Embedded with advanced features
- * Reliability
- * Ergonomics
- * Cost reduction.

This document describes how to manage the DICOM Imager (NANORIP with an PRINTRONIX printer).

1.2 **Document term**

The DICOM conformance statement, license, control panel and web interface messages are all provided in English in this document.

2. **Start**

To perform the operations described in this guide, you must have a NANORIP as well and an PRINTRONIX LP654C or LP844C printer.

This section describes the procedure for configuring and using the NANORIP with an PRINTRONIX printer.

NANORIP Server Configuration Guide for the PRINTRONIX Printer Series

Prerequisites:

- A computer connected to the same network as the printer. (See step 1)
- Graphics software such as Adobe **Photoshop**®, **MS Paint**®, **GIMP**®, etc.... (See step 6)
- Charruasoft **TestSCU** software(See step 6)
- Downloadable from this link <https://www.charruasoft.com/products/free/testscu.zip>

Index:

- Step 1 : Setting up the Printer	Page 8
- Step 2 : Access to the DICOM server home page	Page 12
- Step 3 : Description of [SETUP] menu	Page 14
- Step 4 : Description of [CONFIG] menu	Page 21
- Step 5 : Description of [LOGS] menu	Page 23
- Step 6 : Brochure customization	Page 24

Step 1: Setting up the printer.

Let's start by making sure that the PRINTRONIX printer is correctly installed with all its consumables and paper, so that we can continue with the configuration of the NANORIP server.

CAUTION!


You must set the printer sleep mode with a minimum value of 5 minutes (please see printer documentation).



Once the printer is ready, we will now install NANO DICOM RIP.



Secure NANORIP to the printer back with the fasteners provided as shown in the picture below.



Connect the USB power cable to the NANORIP USB type C port, as well as to the USB type A port of the printer as shown in the photo below.


If your printer does not have a USB Type A port, please use a 5V 2A Power Adapter.




Connect the USB printing cable to the NANORIP USB 3.0 type A port, as well as to the USB type B port of the printer as shown in the photo below.



Connect the NANORIP using a network cable (RJ45) to your network socket as shown in the photo below.



Once NANORIP is plugged to the printer and connected to your network, plug in the printer power cable and power on the printer. Enter using the printer's control panel (see the printer documentation), the IP address, subnet mask and the gateway that will have been provided to you by the person in charge of your computer network.



Step 2 : Access to the DICOM server home page.


Once step 1 has been validated, we can access the reprographer's management interface. To do this, from a computer connected to the same network as NANO DICOM RIP, please launch an internet browser (Internet Explorer, Firefox etc.)



CAUTION!

When starting for the first time, NANORIP may take up to 5 minutes to configure and retrieve printer information.

A web page will be displayed:




PS: If you do not have internet access, no page will appear but the procedure is still doable as long as you can access the address bar (see below)

In your web browser [address bar] enter the printer IP address (the one we has set previously):

<http://<<Printer IP>>>

(For our example, our printer's IP address is 192.168.1.34)


Either in our case: <http://192.168.1.34>



CAUTION!


The operation may fail if performed immediately after [Ready to Print] is displayed on the printer's control panel. Please wait 5 to 10 seconds and then start again.

DICOM configuration page will then be displayed:




NANORIP Home Page.

Step 3 : Description of [DICOM Nodes] menu




To open the Setup menu, click on the [DICOM Nodes] button.

This will show us the list of printing rules represented by their AETitle.



You can add new rules by clicking on [Add new Node], edit them by click on their name (in our case, [NANO]), or delete them by click on [X] in front of their name.



Please note [NANO] AET is mandatory for NANORIP work properly and cannot be deleted.

At any time you can go back to home page by click on [HOME]

Click on NANO to edit this printing rule.

NANORIP	
Current AET	NANO
Calibration Mode	Off
LOG Level	Only Error LOG
Finishing Mode	Standard
Media Type	Paper
Media Size	Automatic
Printer Tray	Automatic
Header	
Footer	Copy for the patient
Number of Copy	Client
Color Mode:	Standard
Contrast:	0
Brightness:	0
Cyan Gamma:	0
Magenta Gamma:	0
Yellow Gamma:	0
Black Gamma:	0
Border Density	Client
Empty Image Density	Client
Sharpen Image Filter	Off

SAVE Reset

Copyright © 2022 RB3L Ltd. All rights reserved.

Description of the different options for the printing rules:

A large number of rules can be added to the system, and for each, the following information is configurable:

Current.AET	NANO
-------------	------

[Current AET] : AET is an abbreviation for Application Entity Title. The invocation method is assumed to conform to the DICOM Print SCU protocol.

AETs defined in this rule have a maximum length of 16 characters (a - z, A - Z, 0 - 9).

Note: • Names are case sensitive.

- the AET is only definable and modifiable when creating a rule.
- In error case, the printer will use the default AET [NANO] and print the overlay text: "This AET (with the name of the AET) **does not exist**".

LOG.Level	Only Error LOG ▼
-----------	------------------

[Log Level] : This function is used to display the DICOM dialog logs when printing. This is useful for identifying the source of an error in case of problems. You can manage 3x logs level, [Only Error LOG], and [High Level Msg].

[Only Error LOG] only display error messages.

[High Level Msg] displays all DICOM Tag messages.

Finishing.Mode	Standard ▼
----------------	------------

[Finishing Mode] : This function allows you to choose the printout in page mode [STANDARD] or brochure mode [BROCHURE].

When brochure mode is enabled [BROCHURE] 3 x new submenus appear.

Finishing.Mode	Brochure ▼
...Brochure.format	A4 ▼
...Add.a.Blank.Page	Off ▼
...Image.on.the.Last.brochure.Page	Off ▼

[Brochure format] : This function allows you to choose the printout paper format for your Brochure [A4], [LETTER], [A3] or [TABLOID]

Note: • [A3] and [TABLOID] is only available when using an PRINTRONIX CLP844C printer.

[Add a Blank Page] : This function allows you to add a blank page behind the cover page to paste a CD-Rom or to staple the examination report.

[Image on the Last brochure Page] : This function allows you to use the last brochure page to put images in place of the end cover page.

Media Size	A4 ▼
------------	------

[Media Size] : This function allows you to select the paper size you want to use, the sizes range from [A6] to [A3 NOBI] depending on the printer used

Printer Tray	Automatic ▼
--------------	-------------

[Printer Tray] : This function allows you to choose the paper tray to use for printing [AUTOMATIC] the printer will select the paper tray according to the requested format.

Header	
Footer	Copy for the Patient

[Header] : This function allows you to customize printout header.

Note: • Customization is a single line of text only.

[Footer] : This function allows you to customize the footer of your prints

Note: • Customization is a single line of text only.

Number of Copy	Client ▼
----------------	----------

[Number of Copy] : This function allows you to choose the number of copies. Copies numbers varies from [1] to [5], [CLIENT] means that the copies number come from the modality.

Color Mode	Standard ▼
------------	------------

[Color Mode]: This function sets the filter action level. Choices are [STANDARD] [ENHANCED] [PHOTO] or [MONOCHROME]. By default, [STANDARD] is selected.

[STANDARD], Sets printing color images with low toner usage (CMYK).

[ENHANCED], Sets printing color images in enhanced quality (CMYK).

[PHOTO], Sets printing color images in highest quality (CMYK).

[MONOCHROME], Sets printing color images in monochrome (K toner only).

Contrast	0 ▼
Brightness	0

[Contrast] : This function allows you to modify printout contrast. The choices range from [-10] to [+10]. Negative values reduce contrast to flatten the image, [0] means no change, positive values increase contrast to improve the image.

[Brightness] : This function is used to change printout brightness. The choices range from **[-10]** to **[+10]**. Negative values reduce light to darken the image, **[0]** means no change, positive values increase light to brighten the image.

Cyan.Gamma:	<input type="text" value="0"/>
Magenta.Gamma:	<input type="text" value="0"/>
Yellow.Gamma:	<input type="text" value="0"/>
Black.Gamma:	<input type="text" value="0"/>

[Cyan Gamma] : This function is used to adjust printout cyan gamma channel. The choices range from **[-30]** to **[+30]**.

[Magenta Gamma] : This function is used to adjust printout magenta gamma channel. The choices range from **[-30]** to **[+30]**.

[Yellow Gamma] : This function is used to adjust printout yellow gamma channel. The choices range from **[-30]** to **[+30]**.

[Black Gamma] : This function is used to adjust printout black gamma channel. The choices range from **[-30]** to **[+30]**.


Border.Density	Client ▼
Empty.Image.Density	Client ▼

[Border Density] : This option defines the background images, in black or white. Choices are **[BLACK]** **[WHITE]** or **[CLIENT]**. **[CLIENT]** means that the settings will be managed by the DICOM calling modality


[Empty Image Density] : This option defines the image-free zone in black or white. Choices are **[BLACK]** **[WHITE]** or **[CLIENT]**. **[CLIENT]** means that the settings will be managed by the DICOM calling modality




Once the values have been edited, you can save this configuration by click on **[SAVE]** or return to the previous values by click on **[Reset]**. Click on **[SAVE]** to saves and sends you to the AETs list page.




Click on [HOME] to go back to the home page.



Step 4 : Description of [Server Config] menu



Click on [Server Config] to enter in the config menu.



Server IP Address	192.168.1.37
Server Subnet Mask	255.255.255.0
Server Gateway	192.168.1.1

[Server IP Address] : Displays the NANORIP IP address.

[Server IP Subnet Mask] : Displays the NANORIP subnet mask.

[[Server IP Gateway](#)] : Displays the NANORIP gateway.

Server DICOM Print Port	<input type="text" value="11112"/>
Server DICOM Store Port	<input type="text" value="11114"/>

[[Server DICOM Print Port](#)] This function is used to display the input DICOM PrintSCP port.

[[Server DICOM Store Port](#)] This function is used to display the input DICOM StoreSCP port.

DICOM Store Temporization:	<input type="text" value="3"/>
Brochure Temporization:	<input type="text" value="5"/>

[[DICOM Store Temporization](#)] : This function is used for StoreSCP to set the waiting time to get all images before sending DICOM printing request. The choices are from [1] to [60] seconds.


[[Brochure Temporization](#)] : This function is used for PrintSCP to set the waiting time to get all images before sending Brochure printing request. The choices are from [1] to [60] seconds.

Printer Model	<input type="text" value="C834"/>
Printer Serial	<input type="text" value="AL54007618"/>

[[Printer Model](#)] : This function is used to display the connected printer model.

[[Printer Serial](#)] : This function is used to display the printer serial number.

Click on [[HOME](#)] to go back to the home page.



NANO DICOM RIP

Control Interface

DICOM Nodes

Server Config

DICOM Logs

R B
3 L

Copyright © 2022 RB3L Ltd. All rights reserved.

Step 6 : Brochure customization.

6.1 Introduction

NANORIP can print your exams as an A4, Letter, Tabloid or A3 brochure with customizable cover and end pages. We will see here how to proceed with these pages creation.

Prerequisites :

- A computer connected to the same network as the printer. (See step 1)
- Graphics software such as Adobe **Photoshop**®, **MSPaint**®, **GIMP**®, etc... (See step 6)
- Charruasoft **TestSCU** software(See step 6)
- Downloadable from this link <https://www.charruasoft.com/products/free/testscu.zip>

6.2 Démarrage

From your computer, using your graphics software, create a Jpeg A4 image size for A4 or A3 brochure or Letter for a Letter or Tabloid brochure in order to use it for your brochure cover page. Repeat the same operation for the brochure end page.

CAUTION!

To use the previously created cover and end pages, we need to send them to the NANORIP server using TestSCU.exe set with specific AETs for the cover and end pages. It is not necessary to create these AETs on the NANORIP server.

Once these operations have been performed, and still from your computer, unzip the testscu.zip file and double-click on TestSCU.exe to start it




At the top of PrintSCU software, there is a list of buttons



Click the setup button to configure the software to connect to the printer.

Enter the printer's DICOM settings in this window as shown in the photo below.



Enter the NANORIP IP address in the Remote IP Address box (for our example, we entered 192.168.0.11)

Enter the NANORIP server port in Remote Port Number box (we entered 11112 which is the NANORIP default port)

Important:

The AETs to be used for brochure cover page are:

[A4-COVER] for A4 and A3 paper size or [LETTER-COVER] for Letter and Tabloid paper size.

The AETs to use for brochure end page are:

[A4-BACK] for A4 and A3 paper size or [LETTER-BACK] for Letter and Tabloid paper size..

Enter the proper AET to use in remote AET box (for our example, we want to create the A3 brochure cover page , we entered A4-COVER)

Once you have finished entering all the parameters, you must click [OK] to save and start using TestSCU.

You are now on the main PrintSCU window.



We are going to test the DICOM connection, click on the echo check button

If the parameters you entered are OK, this window will appear with the message <<ECHO OK>>.




We are now going to send our cover image to the NANORIP server



Click on the DICOM Print button

Select the Jpeg file of your cover image then click [Open].



After a few seconds, the image will be loaded into the NANORIP server and the message <<PRINT OK>> will appear in this window.



Restart the operation to install the brochure end page in NANORIP server.

Your printer is now ready to print in brochure mode with your cover page and end page elements.

Dealer's stamp:



71-75, Shelton Street,
Covent Garden, London,
WC2H 9JQ, ENGLAND

info@rb3l.com
<http://www.rb3l.com>